



**DIAGNOSTIC KIT FOR
DETERMINATION OF
FRUCTOSAMINE ACTIVITY**

Kit name	Kit size	Cat. No
DR.FRUCTOSAMINE	1 x 15 ml	FRU01DR
DR.FRUCTOSAMINE	1 x 30 ml	FRU02DR

SPECIMEN

Use fresh patient serum.
Samples are stable for a week at 2-8 °C or for 6 months at - 20 °C.

INTRODUCTION

The fructosamine are formed in blood from glucose present therein. The carbonyl group of the glucose reacts with free protein amino residues causing the formation of Schiff's base. The half life time of the fructosamine is 17-20 days. So fructosamine determination is suitable for a long-term (1-3 weeks) monitoring of sugar metabolism for patients with diabetes, especially with type II diabetes mellitus and also suitable for drug efficacy monitoring.

PROCEDURE

These reagents may be used both for manual assay and in several automatic analysers.
Applications for them are available on request.

METHOD PRINCIPLE (NBT method)

The serum's fructosamine is one kind of macromolecule alkone amines compounds. It can make nitro blue tetrazolium (NBT method) chloride reduced to formazan under the alkalinity condition. The quantity of the formazan is direct proportional to the fructosamine concentration. The color measured at 540 nm (530-550 nm), is directly proportional to the fructosamine concentration.

Manual procedure

wavelength 546 nm
temperature 37°C
cuvette 1 cm

Pipette into the cuvette:

	Calibrator (C)	Test (T)
R1 Fructosamine Reagent	750 µl	750 µl
R2 Fructosamine Reagent	250 µl	250 µl

Mix well & incubate for 5 min. at 37°C, then add

Calibrator	40 µl	-
Serum	-	40 µl

Mix well and after exactly 180 secs read the absorbance A1 of the Test (T) and Calibrator (C) against reagent blank. In next 120 secs repeat absorbance reading A2 and calculate A (A2-A1) for test and calibrator.

REAGENTS

Package

DR.FRUCTOSAMINE

R1-Fructosamine Reagent	1 x 11.25 ml	1 x 22.5 ml
R2-Fructosamine Reagent	1 x 3.75 ml	1 x 7.5 ml
R3-Fructosamine Calibrator	1 vial	1 vial

The calibrator value has mentioned in the vial label.

Calculation

$$\text{Fructosamine activity } \mu\text{mol/L} = \frac{\Delta A (T)}{\Delta A (C)} \times \text{Calibrator concentration}$$

Working reagent preparation and stability

The reagents are ready to use. The reagents when stored at 2-8°C are stable up to expiry date printed on the package. **The assay kit reagents are stable for 1 month after opening and kept at 2-8°C.**

The reagents are stable for 7 days on board the analyser at approximately 2-8 °C.

REFERENCE VALUES

Serum	≤ 286 µmol/L
-------	--------------

It is recommended for each laboratory to establish its own reference ranges for local population.

Concentrations in the test

Carbonate buffer	0.1 mmol/L
Detergent	1.0 %
Preservative	0.05 %
Carbonate	0.1 mmol/L
Nitrotetrazolium blue chloride	0.5 mmol/L

QUALITY CONTROL

To ensure adequate quality control, each run should include assayed normal and abnormal controls. If commercial controls are not available it is recommended that known value samples be aliquoted, frozen and used as controls.

Warnings and notes

- (i) Products for in vitro diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
- (ii) Reagents contains Sodium Azide, Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
- (iii) Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.
- (iv) All specimens used in this test should be considered potentially infectious. Universal precautions, as they apply at your facility, should be used for handling and disposing of materials during and after testing.

PERFORMANCE CHARACTERISTICS

Linearity: up to 1000 µmol/L. If the sample activity exceeds 1000 µmol/L, dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 2000 mg/dl, intralipid up to 2000 mg/dl, bilirubin up to 30 mg/dl, ascorbic acid up to 10 mg/dl do not interfere with the test.

Specificity / Interferences

Repeated measurement using the same serum samples 10 times, the measured values of the coefficient of variation (CV) should be ≤ 10%.

Inter Assay Precision:

Consecutive three batches kit difference between the grant shall be ≤ 10%.

ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 546 nm
- thermostat at 37°C;
- general laboratory equipment;

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Howley JEA, Browning MCK, Fraser CG, Assay of serum fructosamine that minimizes standardization and matrix problems: Use to assess components of biological variation. Clin. Chem 1987; 33: 269-272.

SYSTEM PARAMETERS

Method	: Fixed Time (2-point)
Wavelength	: 546 nm
Zero Setting	: Reagent Blank
Temperature Setting	: 37°C
Incubation Temperature	: 37°C
Incubation Time	: ----
Delay time	: 180 secs
Read time	: 120 secs
No. of Reading	: 2
Interval time	: ----
Sample Volume	: 0.04 ml (40 ul)
Reagent Volume	: 1.0 ml (1000 ul)
Calibrator Concentration	: Refer Calibrator vial
Units	: umol/L
Factor	: ---
Reaction slope	: Increasing
Linearity	: 1000 umol/L



Genuine Biosystem



A Genuine Service For better Tommorrow

**No.18/128, 3rd Floor,
Shanthy nagar 1st Street,
Chrompet, Chennai - 600044, India.
Ph: +91-44-22651845
Email: genuinebiosystem@gmail.com
website: www.gb-group.co.in**