

N BIO - TRIGLYCERIDES (GPO /PAP method)



KIT NAME	KIT SIZE	CAT. NO
N BIO - Triglyceries	2 x 50 ml	MTGL02050M

INTRODUCTION

Triglycerides are built of glycerol molecule sterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

triglycerides + H2O LPL → glycerol + fatty acids

glycerol + ATP GK→ glycerol-3-phosphate + ADP

glycerol-3-phosphate + O2 GPO → dihydroxy-acetone-phosphate + H2O2

H2O2 + 4-AAP + T005 POD→ quinoneimine dye + 2H2O

The colour intensity is proportional to the triglycerides concentration.

KIT CONTENTS

Reagent name	MTGL02050M
R1 - Triglycerides reagent	2 x 50 ml
R2 - standard	1 vial

Refer standard value mentioned in the vial.

WORKING REAGENT PREPARATION AND STABILITY

The reagent supplied is ready for use.

CONCENTRATIONS IN THE TEST

buffer PIPES (pH 7.2)	50 mmol/l
4-aminoantipyrine (4-AAP)	0.4 mmol/l
ATP	1.0 mmol/l
Chromogen	2.0 mmol/l
glycerol kinase (GK)	≥300 U/L
glycerol-3-phosphate oxidase (GPO)	≥1000 U/L
peroxidase (POD)	≥500 U/L
lipoprotein lipase (LPL)	≥2000 U/L
Activators & Stabilizers.	

WARNINGS AND NOTES

- Product for in vitro diagnostic use only.
- If the reagent is deteriorated, the working reagent absorbance > 0.5 (read against distilled water, wavelength λ=505 nm. And also physically we can observe the turbidity. Moreover it leads failure to recover control values within the acceptable range. Avoid contact with skin and mucous membranes.

ADDITIONAL EQUIPMENT

- Automatic analyser or photometer able to read at 505 nm
- Thermostat at 37°C
- General laboratory equipment

SPECIMEN

- Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis.
- Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values. Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 3 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples.

PROCEDURE

These reagents may be used both for manual assay and in several automatic analysers. Applications for them are available on request.

Wavelength	505 nm
Temperature	37°C
Cuvette	1 cm

Pipette into the cuvettes:

	Blank (B)	Standard (S)	Test (T)
R1 Triglycerides Reagent	1000 µl	1000 µl	1000 µl
Bring upto the temperature of determination. Then add			
Distilled water	10 µl		
R2 - standard		10 µl	
Sample			10 µl

Mix well, incubate for 10 min. at 37°C or 20 min. at 25-30°C. Read the absorbance of test A(T) and standard A(S) against blank (B).

CALCULATION

Triglycerides concentration = A (T) / A(S) x standard concentration

From calculated triglycerides concentration value subtract 0.11mmol/l (10 mg/dl), which corresponds to average amount of free glycerol in serum.

REFERENCE VALUES

25 - 160 mg/dl

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

To Ensure adequate quality control, each run should include assayed normal and abnormal controls. If commercial controls are not available it is recommended that known value samples be aliquoted, frozen and used as controls.

For Fully Automated analyzers by using multicalibrators or triglycerides standard the calibration curve can plot and the same should be prepared every 8 weeks or with change of reagent lot number.

PERFORMANCE CHARACTERISTICS

-Sensitivity / Limit of Quantitation: 2.0 mg/dl (0.023 mmol/l).

-Linearity: up to 800 mg/dl (9.06 mmol/l). For higher triglycerides concentrations dilute the sample wit

Specificity / Interferences Haemoglobin up to 3.75 g/dl, bilirubin up to 20 mg/dl and ascorbate up to 62 mg/l do not interfere with the test.

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Jacobs N.J., Van Denmark P.: J. Arch. Biochem. Biophys. 88, 250-255 (1960).
2. Kodischek L.K., Umbreit W.W.: J. Bacteriol. 98, 1063-1068 (1969).
3. Trinder P.: Ann. Clin. Biochem. 6, 24-27 (1969).
4. Schettler G., Nussel E.: Arb. Med. Soz. Med. Prav. Med. 10, 25 (1975).
6. Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 610, (1995).
6. Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, 2209, (1994).
7. Dembinska-Kiec A., Naskalski J.W.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Volumes, 575, (1998).

SYSTEM PARAMETERS

Method	End point
Wavelength	505 nm
Zero Setting	Reagent blank
Temperature Setting	37°C
Incubation Temperature	37°C
Incubation Time	10 mins.
Delay time	----
Read time	----
No. of Reading	----
Interval time	----
Sample Volume	0.01 ml (10 µl)
Reagent Volume	1.0 ml (1000 µl)
Standard Concentration	Refer standard vial
Units	mg/dl
Factor	----
Reaction slope	Increasing
Linearity	800 mg/dl



Genuine Biosystem Private Limited
Plot No.97 & 98, kattabomman street,
Parvathy Nagar Extension,
Old Perungalathur, Chennai - 600063, India.
Ph: +91-44-48681845
Email: genuinebiosystem@gmail.com
website: www.genuinebiosystem.com