

GBLISA-W

A Rapid Microplate ELISA Washer



- ▶ Compatible with U, V or flat-bottom plates
- ▶ Upto 50 Wash protocols programmable
- ▶ Minimal residual volume
- ▶ Unattached absorb head and dispense head
- ▶ Positive displacement pump with low noise
- ▶ Overfill washing and overflow protection
- ▶ Automatic waste liquid warning
- ▶ Bottom washing function
- ▶ 8-way and 12-way wash heads interchangeable

Genuine Biosystem



A Genuine Service For better Tommorrow